


Choice changes everything!
The C-MAC® HD series



The C-MAC® system – We understand airway management!



The C-MAC® HD series

What do you expect from an ideal airway management system? What modifications can we make that would enable you to work even better with our products? How can we increase our focus on you and your patients? - Our answer: The C-MAC® HD series.

Let our C-MAC® HD video laryngoscopes and our C-MAC® HD monitor take you to the next level of high-definition imaging. Our HD visualization forms the basis for even more realistic views of the patient airway while the advanced technology of our endoscopes ensures that the quality of patient data is higher than ever before.

Discover the benefits of having two endoscopes from the C-MAC® family connected at the same time, either in the Picture-in-Picture (PiP) or Side-by-Side (SbS) mode. You can also save your images and videos in the selected mode.

Make a difference and enjoy the advantages offered by a personalized airway setting for daily routines and challenges or for the training of young physicians.


Your expertise can save lives; benefit from ours, namely, the C-MAC® HD series.



The C-MAC[®] HD Monitor: The



C-MACtouch – Interactive FULL HD touch screen that can also be used while wearing OR gloves




C-MACdual – Picture-in-Picture and Side-by-Side modes



C-MACrec – Choice between Click-to-Record, Auto-Record or Time-Machine mode*



C-MACmic – Integrated microphone allows audio comments to be added to your videos




C-MACsolution – Increased diversity of visualization through maximum flexibility of single-use and reusable products



C-MACtimer – Timer for documenting intubation time for the purposes of training or studying


* Thanks to the C-MAC[®] HD monitor, when the Time Machine mode is activated and the recording button is pressed, 60 seconds of a video sequence is automatically saved before actuation.



C-MAC[®] HD Video Laryngoscopes




C-MAC[®] S Video Laryngoscopes



C-MAC[®] VS (Video Stylet)



Heart of the C-MAC® System



**FIVE S (Flexible Intubation
Videoendoscopes
for single use)**



**FIVE (Flexible Intubation
Videoendoscope)**



**Retromolar
Intubation Endoscopes**



Intubation Fiberscopes



The C-MAC® POCKET MONITOR: The End

- 3.5" OTI® (Open To Intubate) display with wide viewing angle (160°)
- Universal C-MAC® system interface, tailored to your airway management solutions from KARL STORZ
- High-contrast images, even in bright daylight
- Real-time documentation for quality assurance
- Exchangeable, reusable lithium-ion battery – operation time approx. 60 min with capacity indicator in the monitor, uninterrupted readiness
- Protected against dust and continuous immersion (IP68)
- Additional standard: RTCA/DO-160G
- Flexible positioning of the monitor for the intubation of difficult-to-access patients



C-MAC® HD Video Laryngoscopes

Standard MACINTOSH and MILLER video laryngoscopes for safe intubation in both direct and indirect laryngoscopy

D-BLADE and D-BLADE Ped. with strong curvature for the difficult airway



C-MAC® S Video Laryngoscopes

Single-use video laryngoscopes for minimizing the risk of cross-contamination



Emergency Tool for Time-critical Situations



C-MAC® VS (Video Stylet)

Deflectable distal tip combines the advantages of both rigid and flexible intubation endoscopes



C-MAC® PM Connecting Cable

The C-MAC® PM connecting cable enables flexible positioning of the monitor under cramped conditions

The C-MAC® HD Video Laryngoscope

For clearer visualization.

- We raise your standards – visualization of intubation in HD (800p)*
- Lighter, better impact-resistance, more robust – with titanium handle
- When seconds count – plug & play with C-MAC® HD
- BlueButton for real-time documentation in HD*
- Added safety thanks to the improved design of the telescope mount
- Can be combined with 3rd and 4th generation C-MAC® monitors
- Protect yourself and others – maintain your distance with C-MAC® HD video laryngoscopes
- Protected against dust and continuous immersion (IP68)

* only in conjunction with Monitor 8404ZXK



MACINTOSH #4



MACINTOSH #3



MACINTOSH #2



MACINTOSH #0



MILLER #2



MILLER #1



MILLER #0



C-MAC® D-BLADE



C-MAC D-BLADE Ped.


The C-MAC® S IMAGER

A solution with single-use video laryngoscope blades

The single-use video laryngoscope blades work in perfect synergy with the reusable C-MAC® S IMAGER and are designed to reach the highest hygiene standards.

- Blade and handle form one continuous piece: Minimizes the risk of cross-contamination
- Exchange of C-MAC® S video laryngoscopes within seconds
- BlueButton: Documentation with innovative multifunctional button and individual color coding
- Straightforward Plug & Play thanks to the universal C-MAC® system interface
- Protected against continuous immersion (IPX8)

C-MAC® S IMAGER



C-MAC® S Pediatric IMAGER




The Gold Standard: The C-MAC® FIVE for the Anticipated Difficult Airway

Fine, fast, focused – visualization in flexible intubation

With the Flexible Intubation Videoendoscopes (FIVE) 6.5 x 65, 5.5 x 65, 4.0 x 65 and 3.0 x 51.5, KARL STORZ once again sets a new direction in airway management. Similar to the C-MAC® HD video laryngoscopes, the FIVEs deliver clear, pixel-free images without a Moiré effect. The flexible intubation videoendoscopes are directly connected to the C-MAC® HD monitor. This allows you to switch to the video laryngoscope very rapidly whenever necessary.

- Compatible with the C-MAC® HD monitor
- Compact design and ergonomic handle shape
- FIVE 3.0 and FIVE 4.0 for position check of Double Lumen Tubes (DLT)
- For airway inspection
- Practical tube fixation thanks to the special tube holder
- High image resolution and 4:3 video format
- Integrated LED light source
- Image and video documentation via buttons on handle



The C-MAC® Video Stylet

Redefine your limits!

The C-MAC® VS can be considered a successor to the familiar retromolar intubation endoscope. Thanks to its sheath and deflectable tip, the C-MAC® VS combines the advantages of both rigid and flexible intubation endoscopes and is therefore of particular benefit in bariatric surgery or for patients with a restricted mouth opening or cervical spine problems. Thanks to its lever, the distal bend can easily be adapted to specific anatomic conditions. The high-resolution CMOS chip ensures clear, pixel-free images in a 4:3 format without a Moiré effect.

- Straightforward Plug & Play thanks to the universal C-MAC® system interface
- BlueButton: Documentation with innovative multifunctional button and individual color coding
- Protected against continuous immersion (IPX8)
- Angulation up to 60° with attached ETT
- Continuous O₂ flow possible via tube adaptor between the tube and the instrument



The C-MAC® FIVE S

Choice changes everything!

Although this endoscope is designed for single use, you can expect the same standards as you would from any other endoscope from KARL STORZ. With respect to sheath stiffness, precise deflection and image quality in particular, we applied the quality standards of our reusable FIVE to the FIVE S.

- Be inspired by our image quality and combine our reusable FIVE with our FIVE S single-use endoscopes
- Discover the tactile feel of the FIVE S with our new handle design and color coding for size identification
- Very good maneuverability, thanks to the special sheath design with additional rigidity, tailored to your indications in airway management
- Perform at a new level – suction within seconds at the touch of a button with a working channel diameter of 2.2 mm
- Let yourself be impressed with the practical tube fixation on the handle and the good gliding properties when withdrawing the endotracheal tube
- Enhance your workflow flexibility and efficiency and benefit from permanent product availability without complex reprocessing cycles thanks to our sterile packaging with click protection



Accessories: Little things that make work easier

We do not stop working when our products are finished – because we want to make life easier for you...

How are flexible fiberscopes and rigid eyepiece endoscopes connected to the monitor?

With the H1 camera head.



What facilitates the introduction of the endotracheal tube?

The C-MAC® GUIDE.




How can you remove foreign bodies from the oropharynx quickly and easily?

With the modified MAGILL forceps.



In addition to the C-MAC® protective cover and the well-organized C-MAC® PM bag, we can also offer other airway management solutions.




What are the applications of the C-MAC® system?

Bedside use

Emergency use



General anesthesia



ICU



Emergency room



Prehospital use



Instrument Set



8404Z XK

C-MAC® HD Monitor for CMOS Endoscopes, screen size 8" with 1920 x 1200 pixel resolution, two camera inputs, a USB and a HDMI port, optimized user interface via touch screen, video and image capturing in real time on SD card, playback of recorded video clips and still images, data transfer from SD card to USB flash drive possible, splash-proof according to IP54, suitable for wipe disinfection, shock-resistant ABS plastic housing, intelligent power management with rechargeable Li-ion batteries, VESA 75 mounting option, power adaptor for EU, UK, USA and Australia, power supply 110-240 VAC, 50/60 Hz, for use with compatible CMOS videoendoscopes including:

C-MAC® HD Monitor
SD Card Ultra, 16 GB
Protection Cap
VESA 75 Quick Clip



8403X

C-MAC® Connecting Cable, with C-MAC® system interface, length 200 cm, for use with C-MAC® Monitor 8403Z XK/8404Z XK, C-HUB® II 20290301 and C-MAC® video laryngoscope series 8403/8404



8403XA

C-MAC® Extension Cable, for the transmission of analog and digital signals, dustproof according to IP50 standards (not waterproof), length 200 cm, for use with C-MAC® Monitor 8403Z XK/8404Z XK and C-HUB® II 20290301 in conjunction with all endoscopes of the KARL STORZ Office Line (8-pin instruments)



8403XDK

C-MAC® POCKET MONITOR Set, unit with LCD monitor and power supply for all C-MAC® video laryngoscopes, with C-MAC® system interface, screen size 3.5", documentation of images and video sequences saved on internal memory, monitor movable via two rotation axes, rechargeable Li-Ion battery, 1 h operation time, exchangeable battery pack, 2 h charging time, power management with capacity indicator, protection class IPX8, for use with C-MAC® video laryngoscope series 8403/8404 including:

Battery, reusable
USB Data Cable
Protection Cap, for reprocessing



8403XDP **C-MAC® PM Connecting Cable**, to maximize the distance between the C-MAC® POCKET MONITOR and C-MAC® video laryngoscopes, with C-MAC® system interface, length 50 cm, for use with C-MAC® video laryngoscope series 8403/8404, C-MAC® POCKET MONITOR 8403XDK



8403XDL **Charging Unit**, for one rechargeable Battery 8403XDA for C-MAC® POCKET MONITOR 8403XD, with power supply and mains adaptor for EU, UK, USA and AUS, power supply 100-240 VAC, 50/60 Hz, suitable for wipe disinfection, for use with Battery for C-MAC® PM



8403XDA **Battery**, rechargeable Li-ion battery for power supply of C-MAC® PM 8403XD, charging via Charging Unit 8403XDL, optimal operation time of 60 min, suitable for wipe disinfection



8403XDD **USB Data Cable**, USB 2.0 port, for data transfer from C-MAC® PM 8403XD to a computer, length 200 cm



TC100EN* **TELECAM C3**, camera control unit with 2 camera inputs (X-LINE and C-LINE), for use with flexible videoendoscopes and one-chip camera heads (up to FULL HD), with digital Image Processing Module and USB storage option, power supply 100-120 VAC/200-240 VAC, 50/60 Hz, including:
Mains Cord, length 300 cm
DVI-D Connecting Cable, length 300 cm
USB Flash Drive, 32 GB
USB Silicone Keyboard, with touchpad, US



809120 **MAGILL Forceps**, for children, modified, length 20 cm, for the endoscopic removal of foreign bodies, for use with video laryngoscopes of sizes 1 and 2



809125 **MAGILL Forceps**, modified, length 25 cm, for the endoscopic removal of foreign bodies, for use with video laryngoscopes of sizes 2-4



8404BXC **C-MAC® HD Video Laryngoscope MAC #4**, CMOS technology, with MACINTOSH laryngoscope blade, size 4, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403ZXX/8404ZXX and C-MAC® PM 8403XDK

8404BX **Same**, with catheter introduction sizes 16-18 Fr.



8404AXC **C-MAC® HD Video Laryngoscope MAC #3**, CMOS technology, with MACINTOSH laryngoscope blade, size 3, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403ZXX/8404ZXX and C-MAC® PM 8403XDK

8404AX **Same**, with catheter introduction sizes 14-16 Fr.



8404KXC **C-MAC® HD Video Laryngoscope MAC #2**, CMOS technology, with MACINTOSH laryngoscope blade, size 2, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403ZXX/8404ZXX and C-MAC® PM 8403XDK



8404EXC **C-MAC® HD Video Laryngoscope MAC #0**, CMOS technology, with MACINTOSH laryngoscope blade, size 0, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403Z XK/8404Z XK and C-MAC® PM 8403XDK



8404MXC **C-MAC® HD Video Laryngoscope MIL #2**, CMOS technology, with MILLER laryngoscope blade, size 2, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403Z XK/8404Z XK and C-MAC® PM 8403XDK



8404GXC **C-MAC® HD Video Laryngoscope MIL #1**, CMOS technology, with MILLER laryngoscope blade, size 1, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403Z XK/8404Z XK and C-MAC® PM 8403XDK



8404DXC **C-MAC® HD Video Laryngoscope MIL #0**, CMOS technology, with MILLER laryngoscope blade, size 0, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403Z XK/8404Z XK and C-MAC® PM 8403XDK



8404HX **C-MAC® HD Video Laryngoscope D-BLADE, for adults**, CMOS technology, for difficult intubation, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, with catheter introduction sizes 16-18 Fr. for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403Z XK/8404Z XK and C-MAC® PM 8403XDK



8404HXP **C-MAC® HD Video Laryngoscope D-BLADE Ped., pediatric**, CMOS technology, for difficult intubation in children, with C-MAC® system interface, documentation of images and video sequences via the BlueButton, for use with C-MAC® Connecting Cable 8403X for C-MAC® Monitor 8403Z XK/8404Z XK and C-MAC® PM 8403XDK



8403XSI

C-MAC® S IMAGER, with C-MAC® system interface, suitable for manual and machine disinfection up to 65 °C and High-Level Disinfection (HLD) acc. to US standards, for use with C-MAC® Connecting Cable 8403X, C-MAC® Monitor 8403ZXK/8404ZXK, C-MAC® PM 8403XDK and single-use C-MAC® S Video Laryngoscopes 051213-10, 051214-10 and 051216-10



051214-10

C-MAC® S Video Laryngoscope MAC #4, with MACINTOSH laryngoscope blade, size 4, for single use, package of 10x 051214-01, for use with C-MAC® S Imager 8403XS or 8403XSI



051213-10

C-MAC® S Video Laryngoscope MAC #3, with MACINTOSH laryngoscope blade, size 3, for single use, package of 10x 051213-01, for use with C-MAC® S Imager 8403XS or 8403XSI



051216-10

C-MAC® S Video Laryngoscope D-BLADE, laryngoscope blade for difficult intubation, for single use, package of 10x 051216-01, for use with C-MAC® S Imager 8403XS OR 8403XSI



8401DS

C-MAC® GUIDE, guide rod made of stainless steel with atraumatic tip, distal region is adapted to the blade shape of the C-MAC® video laryngoscope D-BLADE, fixation of endotracheal tube with the integrated tube holder possible, package of 10, for use with C-MAC® video laryngoscopes



8403XSPK **C-MAC® S Pediatric IMAGER, Set, with C-MAC® system interface, suitable for manual and machine disinfection up to 65 °C and High-Level Disinfection (HLD) acc. to US standards, including Pressure Compensation Cap 11025E and Leakage Tester 13242XL, for use with C-MAC® Connecting Cable 8403X, C-MAC® Monitor 8403Z XK/8404Z XK, C-MAC® PM 8403XD and single-use C-MAC® S Video Laryngoscopes 051110-10, 051111-10**



051111-10 **C-MAC® S Video Laryngoscope MIL #1,** for single use, package of 10, for use with C-MAC® S Pediatric IMAGER 8403XSP



051110-10 **C-MAC® S Video Laryngoscope MIL #0,** for single use, package of 10, for use with C-MAC® S Pediatric IMAGER 8403XSP



10331B XK **C-MAC® VS, Set, rigid intubation videoendoscope, with deflectable tip, CMOS technology, for use with Connecting Cable 8403X for C-MAC® Monitor 8403Z XK/8404Z XK and C-MAC® PM 8403XDK**

Deflection up/down: 60°/0°
Direction of view: 0°
Angle of view: 100°
Working length: 41 cm
Total length: 60 cm
Distal tip outer diameter: 5.5 mm



10331B2K **Retromolar Intubation Endoscope, Set, with movable eyepiece, for ETT >5.5 mm, 35,000 pixels, autoclavable, with Tube Holder 10331BA, for use with H1 camera head**

Deflection up/down: 40°
Direction of view: 0°
Angle of view: 100°
Working length: 40 cm
Total length: 52 cm
Distal tip outer diameter: 5.0 mm



11304BCXK Flexible Intubation Videoendoscope Set 6.5 x 65,
CMOS technology, with suction valve, for use with C-MAC® Monitor 8403ZXX/8404ZXX

Deflection up/down:	180°/140°
Direction of view:	0°
Angle of view:	100°
Working length:	65 cm
Total length:	94 cm
Working channel inner diameter:	3.0 mm
Distal tip outer diameter:	6.3 mm



11303BNXK Flexible Intubation Videoendoscope Set 5.5 x 65,
CMOS technology, with suction valve, for use with C-MAC® Monitor 8403ZXX/8404ZXX

Deflection up/down:	140°/140°
Direction of view:	0°
Angle of view:	100°
Working length:	65 cm
Total length:	94 cm
Working channel inner diameter:	2.1 mm
Distal tip outer diameter:	5.5 mm



11302BDXK Flexible Intubation Videoendoscope 4.0 x 65, Set,
CMOS technology, with suction valve, for use with C-MAC® Monitor 8403ZXX/8404ZXX


Deflection up/down:	140°/140°
Direction of view:	0°
Angle of view:	100°
Working length:	65 cm
Total length:	93 cm
Working channel inner diameter:	1.5 mm
Distal tip outer diameter:	4.0 mm




11301ABXK Flexible Intubation Videoendoscope Set 3.0 x 51.5,
CMOS technology, for use with C-MAC® Monitor 8403ZXX/8404ZXX

Deflection up/down:	140°/140°
Direction of view:	0°
Angle of view:	100°
Working length:	51.5 cm
Total length:	72 cm
Distal tip outer diameter:	2.85 mm




- 0915612-06 **Flexible Intubation Videoendoscope 5.3 x 65,** STERILE 
- sterile, for single use, package of 6, for use with C-MAC® Monitor 8403ZX in conjunction with E-BOX (TP012*) or C-MAC® Monitor 8404ZX, E-BOX for X-LINK (TC028) and TELE PACK+ (TP101)
- | | |
|----------------------------|-----------|
| Deflection up/down: | 180°/180° |
| Direction of view: | 0° |
| Angle of view: | 90° |
| Working length: | 65 cm |
| Working channel diameter: | 2.2 mm |
| Distal tip outer diameter: | 5.3 mm |



- 091361-06 **Flexible Intubation Videoendoscope 3.5 x 65,** STERILE 
- sterile, for single use, package of 6, for use with C-MAC® Monitor 8403ZXK (with E-BOX TP012*) and C-MAC® Monitor 8404ZXK
- | | |
|----------------------------|-----------|
| Deflection up/down: | 180°/180° |
| Direction of view: | 0° |
| Angle of view: | 90° |
| Working length: | 65 cm |
| Working channel diameter: | 1.2 mm |
| Distal tip outer diameter: | 3.5 mm |

* Product is not available for order.
New customers require C-MAC® MONITOR 8404ZX.



- 0916612-06 **Flexible Intubation Videoendoscope 6.5 x 65,** STERILE 
- sterile, for single use, package of 6, for use with C-MAC® Monitor 8403ZX in conjunction with E-BOX (TP012*) or C-MAC® Monitor 8404ZX, E-BOX for X-LINK (TC028) and TELE PACK+ (TP101)
- | | |
|---------------------------|-----------|
| Deflection up/down: | 180°/140° |
| Direction of view: | 0° |
| Angle of view: | 90° |
| Working length: | 65 cm |
| Working channel diameter: | 3.0 mm |
| Outer diameter: | 6.5 mm |



11301BNK1 Intubation Fiberscope 5.2 x 65, eyepiece version, with suction valve, for use with LED Battery Light Sources 11301D1, 11301D3, 11301DE, 11301DF, TELECAM (PAL) 20212030, TELECAM (NTSC) 20212130 and fiber optic light cables

Deflection up/down:	140°/140°
Direction of view:	0°
Angle of view:	120°
Working length:	65 cm
Total length:	93 cm
Working channel inner diameter:	2.3 mm
Distal tip outer diameter:	5.2 mm



11302BDK2 Intubation Fiberscope 3.7 x 65, eyepiece version, with suction valve, for use with LED Battery Light Sources 11301D1, 11301D3, 11301DE, 11301DF, TELECAM (PAL) 20212030, TELECAM (NTSC) 20212130 and fiber optic light cables

Deflection up/down:	140°/140°
Direction of view:	0°
Angle of view:	120°
Working length:	65 cm
Total length:	93 cm
Working channel inner diameter:	1.5 mm
Distal tip outer diameter:	3.7 mm



11301AAK1 Intubation Fiberscope 2.8 x 65, eyepiece version, with suction valve, for use with LED Battery Light Sources 11301D1, 11301D3, 11301DE, 11301DF, TELECAM (PAL) 20212030, TELECAM (NTSC) 20212130 and fiber optic light cables

Deflection up/down:	140°/140°
Direction of view:	0°
Angle of view:	90°
Working length:	65 cm
Total length:	98 cm
Working channel inner diameter:	1.2 mm
Distal tip outer diameter:	2.8 mm



8404YD **Protective Bag, blue**, for the C-MAC® system, made of water-resistant and sturdy material, washable, bag contains separate compartments for the monitor and three C-MAC® video laryngoscopes, for use with C-MAC® Monitor 8404ZXK, C-MAC® video laryngoscope series 8403/8404



8403YE **Bag for Intubation Set -C22-**, ULM model, made of water-resistant and sturdy material, washable, bag contains two compartments with several holding possibilities for C-MAC® video laryngoscopes with C-MAC® PM and for conventional laryngoscopes, for use with C-MAC® PM 8403XDK, C-MAC® video laryngoscope series 8403/8404 and conventional laryngoscopes

8403YES **Bag for C-MAC® PM and IMAGER**, made of water-resistant and sturdy material, washable, bag with holding possibilities for C-MAC® PM and C-MAC® S IMAGER as well as the C-MAC® S Pediatric IMAGER, for use with C-MAC® PM 8403XDK, C-MAC® S Imager 8403XSI and C-MAC® S Pediatric IMAGER 8403XSPK



11301D1 **Battery Light Source LED for Endoscopes**, with fine screw thread, boost mode for temporary increase in brightness, burning time > 120 min, weight approx. 70 g (without battery), for use with KARL STORZ endoscopes



11301D3 **Battery Light Source LED for Endoscopes**, with coarse thread, boost mode for temporary increase in brightness, burning time > 120 min, weight approx. 70 g (without battery), for use with KARL STORZ endoscopes



11301DE **Battery Light Source LED for Endoscopes**, rechargeable, with click connection, boost mode for temporary increase in brightness, color temperature 5500 K, lithium-ion batteries, charging time 60 min, burning time at 100% brightness 40 min, weight approx. 110 g, suitable for wipe disinfection

11301DF **Same**, with fast screw thread



094129 **Battery Charger Li-Ion**, power supply 100–240 VAC, 50/60 Hz, for charging Rechargeable Battery Box 094124 or Battery Light Source 11301DE/DF, for use with Mains Cord 094127



094127 **Mains Cord**, length 150 cm, for Battery Charger 094129



UG120 **Equipment Cart**, narrow, high, rides on 4 antistatic dual wheels equipped with locking brakes, mains switch on cover, energy beam with integrated electrical subdistributors with 12 sockets, grounding plugs, dimensions in mm (w x h x d):
Equipment Cart: 660 x 1474 x 730 mm
Shelf: 450 x 25 x 510
Caster diameter: 150 mm
including:
Base Module, equipment cart, narrow
Cover, equipment cart, narrow
Beam Package, equipment cart, high
3x **Shelf**, narrow
Drawer Unit with Lock, narrow
2x **Equipment Rail**, long
Camera Holder
2x **Mains Cord**, length 100 cm



UG110 **Equipment Cart**, narrow, small, rides on 4 antistatic dual wheels equipped with locking brakes, mains switch on cover, energy beam with integrated electrical subdistributors with 6 sockets, grounding plugs, dimensions in mm (w x h x d):
Equipment Cart: 660 x 1265 x 730 mm
Shelf: 450 x 25 x 510
Caster diameter: 150 mm
including:
Base Module, equipment cart, narrow
Cover, equipment cart, narrow
Beam Package, equipment cart, small
2x **Shelf**, narrow
Drawer Unit with Lock, narrow
2x **Equipment Rail**, long
2x **Mains Cord**, length 100 cm



8401YB **Clamp**, VESA 75 standard, for fixation of C-MAC® monitor to round profile with diameter 20–43 mm and square profile with diameter 16–27 mm, for use with Monitors 8401ZX/8402ZX/8403ZX/8404ZX



ET11-367154 **Round Tube**, 70 cm x diameter 25 mm, for fixation of C-MAC® Monitor and accessories, for use with 8404YA and COR carts



ET11-367224 **Round Tube**, 50 cm x diameter 25 mm, for fixation of C-MAC® Monitor and accessories, for use with 8404YA and COR carts



UG630 **Bracket**, for mounting standard bars (diameter 25 mm) to equipment rails, package of 2



UG623 **Multifunctional Holder**, 2-part, for mounting poles, diameter 25 mm, package of 2, for use with Equipment Carts UGxxx



8404HS **Holder**, for flexible single-use intubation videoendoscopes, two-part, for use with FIVE S



8404VEP **Videoendoscope Protector**, for mounting to equipment rails, incl. installation accessories



UG914

COR Airway Mobile Stand, high, rides on 4 antistatic dual wheels, with locking brakes, high beam module, incl. cable manager and handle, including:

- Subrack**, for COR mobile stand, high
- Top Cover Monitor Holder**, for COR mobile stand
- Equipment Rail**, for COR mobile stand
- Stainless Steel Tube**, 400 mm
- Monitor Module**, for COR mobile stand
- Drawer**, wide, for COR mobile stand



UG915

COR Mobile Stand for C-MAC®, high, rides on 4 antistatic dual wheels, with locking brakes, high beam module, incl. cable manager, handle, including:

- Subrack**, for COR mobile stand, high
- Top Cover Monitor Holder**, for COR mobile stand
- Equipment Rail**, for COR mobile stand
- Drawer**, wide, for COR mobile stand
- Stainless Steel Tube**, 400 mm
- Cross Tube Adaptor**
- Stainless Steel Round Pipe**, 250 mm




8404YA* **Rolling Stand**, for C-MAC® Monitor, height 1340 cm, with five antistatic castors, cross tube connection and crossbar 250 mm, diameter 25 mm, for positioning the monitor, storage tray dimensions (w x d x h): 300 x 200 x 150 mm, with blade holder and cable winder, for use with ET11-367224 Crossbar 500 mm and ET11-367154 Crossbar 700 mm

* Please note: The instruments displayed are not delivered with the stand.



8404AKH **Wire Tray**, including blade holder for C-MAC® and cable winder, for use with standard rails 10 x 25 mm, including:
Wire Tray
Blade Holder for C-MAC®
Cable Winder



39501LC3 **Wire Tray**, for the cleaning, sterilization and storage of three C-MAC® video laryngoscope blades, with silicone holders and lid, external dimensions (w x d x h): 480 x 250 x 105 mm




8404IFH **Holder for Flexible Endoscopes**, for mounting on round profiles, including installation accessories



11301BC **ProShield Protection Tube**, for flexible telescopes, unsterile, for single use, package of 10, distal closed, for use with Holder for Flexible Endoscopes 29005IFH

Notes



Not all products listed in this document are certified according to Directive 93/42/EEC or Regulation 2017/745/EU on medical devices. For this reason, some products that require certification according to this directive/regulation may not be available in these countries.

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: info@ksea.com

STORZ
KARL STORZ — ENDOSKOPE

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com