


OLYMPUS

Your Vision, Our Future

LARYNGOSCOPE SET

A complete set for larynx surgery.




The complete set of rigid telescope for endoscopic larynx surgery: See more, treat better with NBI.

- Brilliant image quality
- Better definition of margins with NBI and HDTV
- High depth of field – no refocussing required
- New approach for better diagnosis



White light



NBI

Telescopes



4 mm telescopes

	Outer diameter	Direction of view	Field of view	Working length	
A22000 A	4.0 mm	0°	60°	28 cm	
A22001 A	4.0 mm	12°	55°	28 cm	
A22002 A	4.0 mm	30°	57°	28 cm	
A22003 A	4.0 mm	70°	90°	28 cm	
A22004 A	4.0 mm	110°	55°	28 cm	
A22005 A	4.0 mm	30°	90°	28 cm	



5 mm telescopes

	Outer diameter	Direction of view	Field of view	Working length	
WA50372 B	5.4 mm	0°	70°	30 cm	
WA50373 B	5.4 mm	30°	70°	30 cm	
WA50374 B	5.4 mm	45°	70°	30 cm	

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.