

No: 10340-000

ESM 2 Suction module



- The type of activation and the degree of suction of the ESM 2 are adjusted via the interface of the ERBEJET 2
- Activation occurs either simultaneous to the application of the waterjet
- The suction pump can be programmed to continue running for a certain time after waterjet is no longer activated.
- Good visibility of the target tissue

Technical Data

| | |
|------------------|-------------|
| Mains connection | |
| System voltage | 120 - 240 V |
| Power frequency | 50 / 60 Hz |

| | |
|--------------------------------|--------------------------------|
| Additional product information | |
| Communication Bus | {ERBE Communication Bus (ECB)} |
| Packaging unit | 1 Piece |

| | |
|---|-----|
| Standards | |
| Protection category acc. to EN 60 601-1 | I |
| Classification acc. to 93/42 EWG | IIa |
| Type acc. to EN 60 601-1 | CF |