

C-MAC HD Monitor: Endoscopy at Your Fingertips

Getting started in endoscopy has never been so easy



Makes veterinary endoscopy fast, easy & affordable

- Compatible with video otoscopes, rigid telescopes, fiberscopes, single-use flexible endoscopes and selected CMOS videoendoscopes*
- Responsive 8" HD LCD touch screen
- Captures still images and videos direct to SD or USB flash drive
- Works with standard electrical power or with the built-in battery
- Intuitive user interface is simple and easy to use

*Not compatible with GI videoendoscopes

STORZ
KARL STORZ – ENDOSKOPE

C-MAC HD Monitor: Endoscopy at Your Fingertips

The C-MAC™ HD monitor is compatible with a wide range of endoscopes, making it a great solution for various endoscopic applications in small animals. The compact and lightweight design is ideal for use in small spaces, while the integrated rechargeable battery provides greater flexibility for use in the field. Recorded images and videos can be reviewed directly on the monitor, allowing for easy review of findings and for client education. The C-MAC™ HD monitor can also be connected to a secondary monitor via the HDMI port for larger display viewing or for presentation to support staff or clients.

The C-MAC™ HD monitor is an excellent entry level solution for getting started in endoscopy or a great secondary diagnostic system for busy practices.

The C-MAC™ HD monitor is recommended for use with:

Video otoscopes

An ideal solution for exam room video otoscopy to increase client satisfaction, client compliance and practice revenue.

Multi-purpose rigid telescopes

A low-cost solution for getting started with a wide range of endoscopic applications in small & exotic animals.



Single-use flexible scopes

An excellent entry point into flexible endoscopy, including upper airway and urogenital endoscopy in small & exotic animals and sinus endoscopy in horses.

Equine dental telescope

The portable design makes equine dental exams, cleaning and teeth extractions in the field fast and easy.

Ordering information

8404Z XK	C-MAC HD Monitor for CMOS Endoscopes, Set
Optional rolling stand	
8404YA	Rolling Stand, for C-MAC monitor
8404AK	Wire Tray, 30 x 20 x 15 cm